



THE QUEEN'S AWARDS
FOR ENTERPRISE
INNOVATION 2013

The Hydrant Tube Director

A simply brilliant way to set up your Hydrant system so that it is there where you need it at all times.

How to use:

1. Push the tube through the tube director sleeve
2. Use the clips and / or Velcro straps to secure the director to your bed, chair, wheelchair etc.
3. Bend the director as needed so as to position the bite valve where you need it to stay



See our website for full details and pricing

W: www.alphamedicalsolutions.com.au

E: Info@alphamedicalsolutions.com.au

P: 1300 783 747 - ABN 17 434 257 512

PO Box 745 - St Ives - NSW - 2075