

EVAC+CHAIR®

The World's No.1 Stairway Evacuation Chair

ARE YOU PREPARED FOR SAFE EMERGENCY EVACUATION?

The Evac+Chair is the **World's No.1**
Emergency Stairway Evacuation Chair



Model

EVAC+CHAIR 300H

DESCRIPTION

The standard model Evac+Chair 300H features a safe working load of 150kg. The chair has a blue textured finish with a contrasting yellow hammock seat

FEATURES

- Easy to use
- One person operation
- Lightweight
- Made in the UK

INCLUDES

- Wall Hooks
- Photoluminescent Signage
- Dust Cover
- User Guide

SPECIFICATION

Height	1050mm (folded)
Width	520mm (folded)
Depth	210mm (folded)
Weight	9.5kg
Safe Working Load	150kg



ACCESSORIES

RECEPTION NOTICE



USER GUIDE



DUST COVER



COMFY SEAT



TRAINING & SERVICING

KEY TRAINER MASTERCLASS

A 3 hour in-house course for up to 3 members of staff that “trains the trainers” in the safe use of the Evac+Chair and provides guidelines and materials that enable successful trainees to undertake in-house training for other colleagues within their organisation.

OPERATOR TRAINING

A 3 hour in-house training course provides training for up to 5 members of staff in the safe use and operation of the Evac+Chair. On course completion all delegates meeting the training standards are certified and deemed competent for their roles.

EVAC+CHAIR SERVICING

An Evac+Chair should be regularly serviced and maintained to ensure its safe operation. To help clients meet this requirement Safety and Mobility can offer a service contract

SERVICING MAINTENANCE INCLUDES

- + VISUAL INSPECTION
- + REPLACEMENT OF ANY WORKING PARTS
- + TESTING
- + RECORDING AND REPORTING

EVAC+CHAIR®

The World's No.1 Stairway Evacuation Chair



EVAC+CHAIR® is a registered trademark of
Stairway Evacuation Chairs Ltd.

© 2023 Stairway Evacuation Chairs Ltd.

0000000000

EVAC+CHAIR® is a registered trademark of
Stairway Evacuation Chairs Ltd.

© 2023 Stairway Evacuation Chairs Ltd.

0000000000

EVAC+CHAIR® is a registered trademark of
Stairway Evacuation Chairs Ltd.

© 2023 Stairway Evacuation Chairs Ltd.

0000000000